

Medical Alert & VoIP

THE 60-SECOND VERSION

- 1 Older medical alert buttons dial over a copper landline and may fail on an internet phone.
- 2 A pressed button that doesn't reach help is silent — you won't know unless you test.
- 3 A cellular medical alert doesn't use your phone line, so switching phones won't break it.
- 4 Call the monitoring company BEFORE you switch, and confirm how your device connects.
- 5 TEST the button after any change — most companies have a test mode for this.

Stuck? Call or text Bill — real help, no rushing.

Call or text: (330) 200-8042

Book a visit: billshometech.com/book · \$99 flat first visit



Does My Medical Alert Still Work on VoIP? Checking Before You Switch

Will your alert button still reach help after you change your phone?

————— **An Extended Guide · No. 1723** —————

Tech Made Simple. Service Made Personal.

Bill Davis — Bill's Home Tech

(330) 200-8042 · bill@billshometech.com · Atwater, Ohio

Book a visit — billshometech.com/book

What's Inside

Everything in this guide, with page numbers. Skip to whichever chapter answers your question right now.

The 60-Second Version	5
1. Why This Guide Matters So Much	5
2. How a Medical Alert Works	6
3. Why the Phone-Line Kind Can Fail	6
4. Why Cellular Alerts Are Safer Now	7
5. Call the Monitoring Company First	8
6. Always Test After a Change	8
7. How to Tell Which Kind You Have	9
8. The Brands You'll Come Across	9
9. What to Ask When Shopping	10
10. A Word About Fall Detection	11
11. Signal and Placement	11
12. Mobile Alerts for an Active Life	12
13. If You're Keeping a Phone-Line Alert	12
14. The Outage Question	13
15. Medical Alert and 911 Together	13
16. Helping a Parent From a Distance	14
17. Testing Made a Habit	14
18. Common Mistakes I See	15
19. A Few Plain Words to Know	15
20. Questions I Get All the Time	16

21. Why the Industry Moved to Cellular	16
22. Don't Wait for the Switch to Check	17
23. What the Monitoring Fee Covers	17
24. Wearing It — Even in the Shower	18
25. Keeping the Pendant Charged	19
26. What Happens When You Press It	19
27. Setting Up a New Cellular Unit	20
28. If You Live Alone	20
29. When to Just Call Me	21
30. What I'd Tell My Own Neighbor	21
Everyday Cheat Sheet	22
A Picture to Follow Along	23
A Second Look	24
A Note From Bill	24
Your Visit — and My Promise	25

SUMMARY

The 60-Second Version

In a hurry? Here's the whole thing on one page. A lot of older medical alert systems were built to call the monitoring center over a copper landline. When you switch your home phone to the internet — cable phone or VoIP — that older button can fail quietly: you press it, but the signal never reaches the center, or arrives garbled. The safest kind of medical alert today is a cellular one that doesn't use your phone line at all. Before you change your phone service, call your medical alert company and ask. And after any change, test the button — the companies have a test mode for exactly this.

- Older medical alert buttons dial over a copper landline and may fail on an internet phone.
- A pressed button that doesn't reach help is silent — you won't know unless you test.
- A cellular medical alert doesn't use your phone line, so switching phones won't break it.
- Call the monitoring company BEFORE you switch, and confirm how your device connects.
- TEST the button after any change — most companies have a test mode for this.

WATCH OUT

This is the one that catches people. A medical alert that worked fine for years can stop reaching help the day you change your phone — with no warning light, no error, nothing. Checking takes one phone call. Please make it.

1. Why This Guide Matters So Much

A button that doesn't reach help is worse than no button.

Of all the things that can go wrong when a home phone changes, this is the one that worries me most, and I want to be straight with you about why. A medical alert button — the pendant or wristband you press to summon help — is only as good as its ability to actually reach the monitoring center. If that connection quietly breaks, you're left with a false sense of safety, which is more dangerous than knowing you have nothing.

Here's the trouble: many older medical alert systems were designed to dial out over a copper landline. As the old copper gets retired and people move to internet-based phone, those older buttons can stop working — and they often do it silently. No alarm, no warning. You'd only find out at the worst possible moment.

So this guide is about one thing: making sure your medical alert still reaches help after any change to your phone. It's factual and practical. I'm not a doctor and I won't give medical

advice — for anything about your health or your device's medical settings, that's a conversation for your doctor and your monitoring company. My job is the technology piece, and I take it seriously.

TIP

If you only do one thing from this guide, call the company that monitors your medical alert and ask: "If I change my home phone, will my button still reach you?" Their answer tells you everything.

2. How a Medical Alert Works

The button, the base, and the monitoring center.

Let me explain how these systems work, because understanding the pieces makes the risk obvious. A typical in-home medical alert has three parts: a button you wear as a pendant or wristband, a base station that sits in your home, and a monitoring center staffed by people around the clock. When you press the button, it signals the base station, and the base station contacts the monitoring center.

At the center, a trained operator comes on the line, talks to you through the base station's speaker, and decides what help you need — a call to a family member, or emergency services to your door. The whole thing hinges on that one link: the base station being able to reach the monitoring center. Break that link, and pressing the button does nothing.

GOOD TO KNOW

The weak point is always the connection between the base station and the monitoring center. How that base station reaches the center — over your phone line or over a cell signal — is the entire safety question in this guide.

3. Why the Phone-Line Kind Can Fail

What happens when copper becomes internet.

The older style of medical alert base station reaches the monitoring center by dialing out over your home phone line, just like a phone would. For decades that meant the copper landline, which carried those signals cleanly. It worked beautifully, and millions of people relied on it.

The problem comes when the home phone moves to the internet. An internet or cable phone carries calls as data, and in the process it can squeeze and reshape the special tones a medical alert uses to talk to the monitoring center. Those tones can arrive garbled, or not arrive at all. The button press seems to go out, but the center never gets a clean signal it can understand.

WATCH OUT

The failure is usually silent. The base station may light up and seem to “call” as always, but the signal reaching the monitoring center is garbled or missing. You cannot tell by looking — you can only tell by testing with the center.

4. Why Cellular Alerts Are Safer Now

They skip the phone line entirely.

Here's the good news, and it's genuinely good. The medical alert industry has largely moved to cellular systems — base stations and pendants that reach the monitoring center over a cell signal, the same kind of network a cell phone uses, with no phone line involved at all. Because they don't touch your home phone, changing or even dropping your home phone doesn't affect them.

A cellular medical alert also brings a safety bonus: many have their own built-in battery, so they can keep working through a power outage, when an internet phone would be dead. For someone who lives alone or has health concerns, that independence from the home's power and phone is exactly what you want in a lifeline.

GOOD TO KNOW

Most medical alert companies now sell cellular systems by default, precisely because the landline is going away. If you're shopping for a new one, a cellular model sidesteps the whole phone-line problem from the start.

5. Call the Monitoring Company First

Before you change a thing about your phone.

If you take nothing else from this guide, take this: before you switch your home phone to anything new, call the company that monitors your medical alert. Tell them exactly what you're planning — “I'm moving my home phone to cable” or “I'm getting an internet phone.” They deal with this every day and can tell you precisely what it means for your device.

- 1 Find the monitoring company's phone number — on the base station, a sticker, or your paperwork.
- 2 Call and tell them you're changing your home phone service.
- 3 Ask: “Will my device still reach you? Do I need a different one?”
- 4 If they recommend a cellular unit, ask how to get one and what it costs.
- 5 Schedule a test for after your phone change is complete.

This one phone call can prevent the worst kind of surprise. The monitoring companies would far rather hear from you ahead of time than have a button fail when it's needed. There's no charge to ask, and it's the responsible thing to do.

TIP

Make this call before you schedule the phone switch, not after. If you need a different device, you'll want it in hand and tested before your old phone line goes away — not scrambling in the gap.

6. Always Test After a Change

The step that proves it actually works.

Asking the monitoring company is step one. Testing is step two, and it's just as important, because the only way to truly know the button reaches help is to try it — safely, through the company's test process. Every reputable medical alert company has a test mode exactly for this, so you can press the button and confirm the signal lands without triggering a real emergency response.

The usual way is to call the monitoring center, tell them you'd like to run a test, and then press your button when they say they're ready. They'll confirm they received a clean signal and that your information came through correctly. Do this right after any change to your phone — and honestly, it's worth doing every month regardless.

WATCH OUT

Never test by pressing the button and waiting to see what happens without telling the company first — that triggers a real response and may send emergency help. Always call the monitoring center and ask for a test so they know it's a drill.

7. How to Tell Which Kind You Have

Phone-line or cellular — and how to check.

You may not know offhand whether your medical alert dials over the phone line or uses a cell signal. That's fine — it's not always obvious, and guessing is exactly what we want to avoid. There are a few ways to find out for sure.

- Look at the base station: a phone-line type usually has a cord running to a telephone wall jack.
- A cellular type typically has no phone cord, just a power cord, and may mention “cellular” or “4G/LTE” on it.
- Check your paperwork or the monthly bill, which often names the type.
- Simplest of all: call the monitoring company and ask them directly.

When in doubt, the phone call to the monitoring company settles it in a minute. They can see exactly what kind of device you have and how it connects to them. Don't rely on a guess for something this important.

GOOD TO KNOW

If your base station has a cord plugged into a telephone wall jack, treat it as a phone-line type until the company tells you otherwise — and that means it needs checking before you change your phone.

8. The Brands You'll Come Across

Named even-handedly — there's no single best.

People often ask me which medical alert brand to get, and my honest answer is that there's no single best — there's the one that fits your needs, your home, and your budget. I don't sell any of them and I don't earn anything from recommending one over another. But it helps to know the names you'll run into so the ads and reviews make sense.

- Life Alert — long-established, widely advertised in-home and mobile systems.
- Medical Guardian — a range of in-home and on-the-go cellular options.

- MobileHelp — known for cellular systems that don't need a landline.
- LifeStation — in-home and mobile alerts with monitoring.
- Bay Alarm Medical — in-home, mobile, and in-car options.

Most of these companies now offer cellular systems as their standard, precisely because the landline is going away. Whichever you consider, the questions are the same: how does it reach the monitoring center, does it have a battery for outages, and what does the monthly monitoring cost. Compare on those, not just the advertising.

GOOD TO KNOW

Prices and features change often, so I won't print numbers that will be stale by next season. Compare each company's current offering, and lean on how the device connects and whether it has outage backup — those matter more than the logo.

9. What to Ask When Shopping

The questions that actually matter.

If you're choosing a new medical alert — whether because your old one won't survive the phone change, or you're getting one for the first time — here are the questions I'd want answered before signing up. They cut through the advertising to what keeps you safe.

- Does it connect over a cell signal, so it doesn't depend on my home phone?
- Does it have a backup battery, and how long does it last in a power outage?
- Is there a monthly monitoring fee, and what exactly does it cover?
- Does it offer fall detection, and is that extra?
- Is there an easy test mode, and how often should I test?
- What is the range of the button around and outside the house?

Write the answers down and compare across a couple of companies. The right choice is the one whose answers fit how you live — whether you're mostly home, out and about, worried about falls, or living alone. There's no prize for the fanciest system, only for the one you'll actually wear and that reliably reaches help.

TIP

The best medical alert is the one you wear every day. A simple pendant you keep on beats a fancy system left in a drawer. Comfort and habit matter as much as features.

10. A Word About Fall Detection

Useful, but understand what it is.

Many medical alerts now offer fall detection — the pendant senses a hard fall and automatically contacts the monitoring center even if you can't press the button. For someone at risk of falls, that can be reassuring, and it's worth asking about. It usually costs a little extra on the monthly fee.

Two honest notes, and I'll stay in my lane here since I'm not a medical professional. Fall detection isn't perfect — it can miss a gentle slide to the floor, and it can occasionally trigger on a hard sit. And it doesn't replace pressing the button when you can. Whether it's right for you is a question for you, your family, and your doctor — not for me. My part is simply making sure whatever device you choose connects reliably.

GOOD TO KNOW

Fall detection is a helpful extra for some people, but it's a supplement to the button, not a replacement. Press the button whenever you can; let the automatic feature be the backstop for when you can't.

11. Signal and Placement

A cellular alert needs a usable signal, like a cell phone.

A cellular medical alert reaches the monitoring center over a cell signal, which means it needs a usable signal where you live — the same as any cell phone. In most of the county that's no trouble, but in a few spots, or inside a thick-walled basement, the signal can be weak. It's worth confirming when the system arrives.

When you set up a cellular alert, test it from the rooms where you spend time, including the bedroom and bathroom, which is where falls often happen. The monitoring company's test mode lets you confirm the signal reaches them from each spot. If one room is weak, the company may have a solution, or a different carrier's system might suit your home better.

WATCH OUT

Test the alert from the bathroom and bedroom, not just the living room. Those are the rooms where help is most often needed, and they're sometimes the ones with the weakest signal.

12. Mobile Alerts for an Active Life

Help that goes where you go.

If you're still out and about — running errands, gardening, walking the neighborhood — there are mobile medical alerts that work anywhere there's a cell signal, not just at home. They use the cell network and often have location built in, so if you press the button at the store or on a walk, the monitoring center can send help to where you actually are.

These are naturally cellular, so they sidestep the phone-line problem entirely. For an active person, a mobile alert can make more sense than a home-only base station. Some people even carry a mobile unit and keep a home base too. As always, how it reaches the center and whether it holds a charge are the questions that matter.

GOOD TO KNOW

A mobile medical alert is worth considering if you leave the house regularly. Because it's cellular by nature, it won't be affected by any change to your home phone — one less thing to worry about.

13. If You're Keeping a Phone-Line Alert

What to do if you can't switch to cellular yet.

Maybe you have a phone-line medical alert you'd rather not replace right now. That can be alright — but only if you've confirmed it still reaches the monitoring center over your new phone, and you keep testing it. Some internet phone setups can carry an alert's signals well enough; others can't. The only way to know is to ask the company and test, not to hope.

If you go this route, treat testing as non-negotiable: run a test with the monitoring company right after the phone change, and then monthly. And keep in mind the outage weakness — a phone-line alert on an internet phone will be dead in a power or internet outage, when a cellular alert with its own battery would still work. Weigh that honestly.

WATCH OUT

A phone-line alert kept on an internet phone must be tested after the switch and regularly thereafter. If the monitoring company can't confirm a clean signal, don't rely on it — move to a cellular unit.

14. The Outage Question

Will the button work when the power's out?

Think through the storm scenario, because that's often exactly when an older person is home alone and most at risk. A phone-line medical alert running on an internet phone goes dead the moment the power or internet does, because the phone it depends on goes dead. That's a serious gap for a device whose whole purpose is emergencies.

A cellular medical alert with its own battery is far sturdier here. It keeps working through a power outage, reaching the monitoring center over the cell network, independent of your home's power and internet. For a lifeline device, that independence is worth a great deal. It's one of the strongest reasons the industry moved to cellular.

WATCH OUT

If the medical alert is important to your household's safety, a cellular unit with a battery is the setup I'd steer you toward — it's the one that keeps working when the power and phone don't. Always keep a charged cell phone as well.

15. Medical Alert and 911 Together

How they fit — and why you still keep a cell phone.

A medical alert and 911 work alongside each other. When you press the button, the monitoring center's operator talks with you and can summon emergency help — effectively reaching 911 on your behalf and staying with you until it arrives. For many people, especially if speaking or dialing is hard, that human on the line is a real comfort.

But a medical alert doesn't replace keeping a charged cell phone. The pendant might be out of reach, out of battery, or out of range when you need it; a cell phone is a second, independent way to reach help. Layering a cellular medical alert and a charged cell phone gives you two separate paths to emergency services — and two is far safer than one.

WATCH OUT

Never treat a medical alert as your only way to reach help. Keep a charged cell phone too. If you ever can't reach either, and you can get to a phone, 911 is always there — on any phone, no plan required, as long as it has a charge and a signal.

16. Helping a Parent From a Distance

What adult children should check.

So many of my calls about this come from adult children worried about a parent living alone. If you're helping from out of town, the medical alert deserves a place on your checklist whenever the parent's phone changes — because they may not realize the button could be affected, and they may not think to tell you.

- Confirm whether the alert is phone-line or cellular — call the monitoring company.
- If the home phone is changing, make sure the alert is checked before, not after.
- Push for a cellular unit with a battery if the current one is phone-line.
- Arrange a monthly test, and ask the company to flag you if it ever fails a test.
- Make sure a charged cell phone is also within the parent's reach.

It's an awkward subject to raise with a parent sometimes, but framing it as “I want to be sure your button reaches help” usually lands gently. Most people would rather know their lifeline works than find out the hard way.

17. Testing Made a Habit

Monthly, and after anything changes.

Testing isn't a one-time thing. Signals drift, batteries age, and phone services get tweaked. The safest households test their medical alert regularly — once a month is the common advice — and always after any change to the phone or internet. It takes a few minutes and it's the only way to be sure.

- 1 Call the monitoring center and tell them you'd like to run a test.
- 2 When they're ready, press your button and let them confirm they received it.
- 3 Confirm they have your correct address and contact information on file.
- 4 Note the date so you remember when you last tested.
- 5 Repeat after any phone change, a move, or a new base station.

TIP

Pick a monthly date you'll remember — the first of the month, say — and make the test part of the routine, like taking out the trash. A tested button is a trusted button.

18. Common Mistakes I See

The slip-ups that put people at risk.

After helping a lot of families with this, the same few mistakes come up again and again. None of them are anyone's fault — the risk is easy to miss — but knowing them ahead of time keeps a loved one safe.

- Assuming the medical alert “just works” on a new phone — it may not.
- Changing the phone service without calling the monitoring company first.
- Never testing, so a silent failure goes unnoticed for months.
- Keeping a phone-line alert that dies with the power, instead of a cellular one.
- Relying on the alert alone, with no charged cell phone as a second path.

WATCH OUT

The most dangerous mistake is the quiet one: a button that looks fine but no longer reaches help. The only cure is to ask the monitoring company and test. Please don't skip it.

19. A Few Plain Words to Know

The jargon, translated.

The medical alert world has its own terms. Here's the short translation so nothing throws you.

- Base station — the box in your home the pendant signals.
- Monitoring center — the staffed center that answers when you press the button.
- Phone-line (landline) unit — an alert that dials out over your home phone line.
- Cellular unit — an alert that reaches the center over a cell signal, no phone line.
- Fall detection — a feature that contacts the center automatically on a hard fall.
- Test mode — a way to press the button and confirm it works without a real response.

TIP

Keep this little list with your guide. When the monitoring company uses one of these words, you'll know exactly what they mean.

20. Questions I Get All the Time

Quick, plain answers.

Will my medical alert stop working when I switch phones? It might, if it's a phone-line type. A cellular one won't be affected. Call the monitoring company to be sure.

How do I test it safely? Call the monitoring center, say you want to run a test, and press the button when they're ready. Never test without telling them first.

Do I need a home phone for a medical alert? No — a cellular alert works with no home phone at all.

Is a cellular alert better? For most people today, yes — it skips the phone line and usually has a battery for outages. But the right one depends on how you live.

GOOD TO KNOW

If your question isn't here, it's probably one I've answered before — just call and ask. For anything about your health or the medical side, though, that's a question for your doctor and the monitoring company.

21. Why the Industry Moved to Cellular

The quiet reason your options look different now.

If you shopped for a medical alert years ago and look again now, you'll notice the choices have shifted — nearly everything is cellular, and the old phone-line units are fading from the catalogs. That's not a coincidence or a sales gimmick. It's the direct result of the copper phone lines being retired across the country. The companies saw that a device depending on a landline was depending on something that's going away, so they moved their whole approach to the cell network.

That shift is genuinely good for you. A cellular alert is simpler to install — no phone jack needed — it keeps working in a power outage thanks to its battery, and it isn't tripped up by the move to internet phone. In other words, the industry already solved the problem this guide is about. The only thing left for you to do is make sure your own device is on the right side of that line, and that's a single phone call to find out.

If your device is older and still ties to the phone line, don't take that as a crisis — take it as a nudge. You're simply due for an upgrade that most people have already made, and it will

leave you with a sturdier lifeline than you had before. That's a good outcome, and an easy one to reach.

GOOD TO KNOW

If a company is still trying to sell you a phone-line medical alert in this day and age, ask why, and ask for the cellular option instead. The cellular unit is where the industry has sensibly landed.

22. Don't Wait for the Switch to Check

Even if your phone isn't changing yet.

Here's a thought worth sitting with. You don't have to be switching your phone today for this to matter. If you have an older phone-line medical alert, it's worth confirming it still reaches help even on your current line — because copper is being retired gradually, and sometimes the underlying service changes without much fanfare. A device that worked when it was installed years ago deserves a fresh check.

So whether or not a phone change is on your calendar, a quick call to the monitoring company and a test is never wasted. Think of it like checking the smoke detectors — you don't wait for a fire to find out whether they work. The same calm, once-in-a-while check keeps your lifeline trustworthy.

GOOD TO KNOW

If it's been a long time since you tested your medical alert, there's no better day than today. One call to the monitoring company tells you whether your button still reaches help — regardless of what your phone is doing.

23. What the Monitoring Fee Covers

Where the monthly money actually goes.

Most medical alerts carry a monthly monitoring fee, and it's worth understanding what you're paying for so you can compare honestly. The fee pays for the staffed monitoring center — real people available around the clock to answer when you press the button, talk you through what's happening, and summon help. That human on the line, any hour of the day, is the heart of what you're buying.

Prices and plans vary a good deal between companies and change often, so I won't print numbers that will be stale by next season. What I'll say is to look at the all-in monthly cost, ask whether the equipment is included or extra, watch for long contracts and cancellation terms, and ask whether features like fall detection or a mobile unit add to the fee. Compare on the whole picture, not the headline rate.

TIP

Ask whether there's a long-term contract or whether you can pay month to month, and what happens if you need to cancel. A good company won't lock you into something you can't leave if your needs change.

24. Wearing It — Even in the Shower

The bathroom is where falls often happen.

A medical alert only helps if you're actually wearing it when something happens, and the hard truth is that many falls happen in the bathroom — getting in or out of the shower or tub, or up from the toilet. So the pendant needs to be something you keep on, and ideally one that's water-resistant enough to wear in the shower.

When you choose a device, ask whether the button is rated to be worn in the shower, because that's exactly where you may need it most. And pick a style you'll genuinely keep on — a comfortable pendant or wristband you forget you're wearing beats a bulky one that ends up on the nightstand. I'll stay out of the health judgments; I just want whatever you wear to reliably reach help.

WATCH OUT

Don't take the pendant off to shower if it's water-resistant — the shower is one of the most common places for a fall. Confirm with the company that your button is safe to wear in the water.

25. Keeping the Pendant Charged

A dead button helps no one.

Different medical alerts handle power differently. Some pendants have a long-life battery the company replaces for you; others, especially mobile units, need charging every so often, like a cell phone. Know which kind you have, because a pendant with a dead battery is just a piece of plastic when you need it.

If yours needs charging, build it into a routine — set it on its charger at the same time each night, or whenever you sit down to a meal. Many systems warn you with a light or a beep when the battery runs low; learn what that warning looks and sounds like so you don't miss it. The base station usually has its own backup battery too, which is what carries it through a power outage.

WATCH OUT

If your pendant needs charging, pair it with a daily habit you already have — charging it while you eat breakfast, say — so it never runs flat without you noticing.

26. What Happens When You Press It

So there are no surprises in the moment.

It helps to know exactly what happens when you press the button, so it's familiar if you ever need it for real. When you press, the base station or mobile unit connects to the monitoring center, and within moments an operator comes on the line and speaks to you through the device's speaker — you don't have to hold anything to your ear.

The operator will ask what's wrong and what you need. Depending on the situation and your instructions on file, they may call a family member or neighbor, or send emergency help straight to your door. They stay on the line with you until help arrives. Because your address and contacts are on file, they know where to send help even if you can't speak clearly — which is exactly why keeping that information current matters.

GOOD TO KNOW

Make sure the monitoring company has your current address and an up-to-date list of who to call. Review it once a year, and whenever a contact's number changes, so help always goes to the right place and the right people are notified.

27. Setting Up a New Cellular Unit

The day-one steps that matter.

If you're moving to a cellular medical alert, the setup is usually straightforward, but a few first-day steps make sure it's truly ready. Don't consider it done until you've worked through them with the monitoring company.

- 1 Plug in the base station and let it find the cell signal; confirm it shows it's connected.
- 2 Confirm the company has your correct address and contact list.
- 3 Charge the pendant fully, and learn its low-battery warning.
- 4 Run a test from several rooms, including the bathroom and bedroom.
- 5 Ask how often to test going forward, and note the date.

Once those are done, the device is genuinely ready, and you can wear it with confidence. If any step is confusing — especially confirming the signal or running the test — that's exactly the sort of thing I can handle in a few minutes on a visit.

TIP

Do the room-by-room test on day one, while someone's there to help. It's much easier to catch a weak-signal room at setup than to discover it in an emergency.

28. If You Live Alone

Building the safety net a little stronger.

If you live by yourself, a medical alert can be one of the most reassuring things in the house — for you and for the family who worry about you. But it works best as part of a small, layered plan, so that no single failure leaves you unable to reach help.

- Use a cellular medical alert with a battery, so it survives a power outage.
- Wear the pendant all day, including in the shower if it's water-resistant.
- Keep a charged cell phone within reach as a second, independent path to help.
- Make sure your home phone's 911 address is registered, as a third path.
- Test the alert monthly, and keep your address and contacts current with the center.

None of this is meant to worry you — it's the opposite. A few small, layered choices mean you never have to wonder whether you can reach help. This is the kind of setup I'm always glad to put together for someone living on their own, and to walk a worried family member through.

29. When to Just Call Me

The technology piece is exactly what I help with.

I wrote this guide so you'd understand the risk and know what to check. The connection side of it — figuring out whether your alert is phone-line or cellular, making sure your new phone won't break it, helping you set up and test a cellular unit, and getting a battery backup in place — is exactly the kind of visit I do. What I won't do is give medical advice; the health decisions stay with you, your family, and your doctor.

I come to you, here in Portage County and the nearby towns. The first in-home visit is a flat ninety-nine dollars, and if a problem can be handled without a trip, I offer remote help too. If you're not completely happy with a visit, I'll make it right or refund it. That's the whole deal — no fine print.

TIP

When you call, my assistant Sophie may answer and take down the details so I can call you back — she's there so you always reach a friendly voice, not a dead end.

30. What I'd Tell My Own Neighbor

The short, honest version.

After 30+ years of professional experience helping folks through exactly this kind of change, here's the plain advice I'd give over the fence. Before you touch your home phone, call the company that monitors your medical alert and ask whether your button will still reach them. If it's an older phone-line type, move to a cellular one that doesn't depend on the phone line at all. And after any change, test it — and keep testing it monthly.

That's really the whole thing. A medical alert is meant to be the one device you can count on in your worst moment. A little checking now makes sure it stays that way — reaching help the instant you press it, power outage or not. For anything about your health, lean on your doctor and the monitoring company; for the technology, lean on me.

GOOD TO KNOW

If any of this feels like a lot, that's normal. You can hand the technology side to me — checking the device, setting up a cellular unit, the battery, the testing — and I'll make sure your lifeline truly reaches help.

CHEAT SHEET

Everyday Cheat Sheet


Keep this by the phone. I want to... Do this. Know if my alert is at risk — call the monitoring company and ask whether it dials over the phone line or uses a cell signal. Protect it before switching phones — call the monitoring company **BEFORE** you change your home phone; move to a cellular unit if yours is a phone-line type. Prove it works — after any change, call the center, ask for a test, and press the button when they're ready; never test without telling them. Stay covered in an outage — choose a cellular alert with a battery, and keep a charged cell phone too. For health questions — call your doctor and the monitoring company. Get tech help — call Bill: (330) 200-8042.

AT A GLANCE

A Picture to Follow Along

Will My Alert Work After I Switch?

Follow the path for your device



Never assume. A pressed button that doesn't reach help is the worst surprise.

After any phone change, test the button — most companies have a test mode.

AT A GLANCE

A Second Look

Two Kinds of Medical Alert

How it connects is the whole safety question

	Phone-line alert	Cellular alert
How it reaches help	Dials your home line	Over a cell signal
Risk on an internet phone	May fail quietly	Not affected
Needs a home phone?	Yes	No
Works in a phone outage?	No	Often, with its battery

Confirm your device's type with the monitoring company; don't guess.

BEFORE YOU GO

A Note From Bill

A quick note from Bill. I write these guides to make home tech a little less frustrating, and I do my best to keep every step accurate and up to date. Still, every home, device, and setup is a bit different, so I can't promise a guide will match yours exactly or fix your particular problem. These guides are for general information only — they're not professional, legal, or financial advice, and they don't replace your device's own manual or the manufacturer's instructions. Please go slow, use your own good judgment, and use everything here at your own risk. If a step involves electrical wiring, a gas appliance, a water heater, plumbing, or climbing a ladder — or if anything just doesn't feel safe — please stop and call a licensed professional (or give me a call at 330-200-8042). Bill's Home Tech and Bill Davis can't be responsible for any damage, data loss, or injury that results.



HOW I CAN HELP

Your Visit — and My Promise

In-home visit — \$99 flat for your first visit. No contracts and no surprise charges, and a 30-day follow-up is included.

Remote Support — \$49 for a screen-share session. I only connect with your permission, you watch the whole time, and you can end anytime. If I can't fix it over the screen-share, I refund the \$49 — and if it turns out you need me in person, that \$49 goes toward your visit. It's included free for Support Plan members (\$39/mo).

GOOD TO KNOW

My promise — If you're not 100% satisfied with your visit, I'll come back at no charge, or refund the full \$99. No questions, no awkwardness. That's how I'd want to be treated.

Call or text **(330) 200-8042** · Book a visit: billshometech.com/book · bill@billshometech.com

This free guide is for general information only — not professional advice, and no substitute for your device's own manual. Any prices or outside products mentioned can change; please confirm current details before you buy. Please use your own good judgment and follow it at your own risk — for anything involving electricity, gas, water, or heights, stop and call a licensed professional. Bill's Home Tech and Bill Davis can't be responsible for any damage, data loss, or injury that results. © 2026 Bill's Home Tech, LLC.